

Collection Protocol DXA_v1

DXA data was collected in the following manner:

Body composition was captured on a GE Lunar DXA or iDXA machine. Participants who had on clothing with snaps, buttons or zippers undressed to their underwear, donned 2 hospital gowns and attempted to void. Those who had light weight clothing without snaps, buttons or zippers were allowed to keep their clothes on. Shoes were removed. Participants were then weighed and laid down on the DXA table according to manufacturer instructions. GE software calculated all tissue masses. In some cases, for large individuals, lateral masses were estimated based on the other side of the body. These estimates are marked with an 'e'.

DXA_v1 was implemented on Aug 15, 2011